

E-Z Floss® Dental Floss Holder

Instructional story to wind "Mr Flossman"™

"Mr Flossman"™ is your friend and wants to help you floss your teeth. But first he needs help from you. "Mr Flossman"™ likes to go skiing, and he needs your help to get him dressed for the cold.

- Using a 12" -15" length of floss
- Hold one end of the floss between his arms (on the X) with your thumb
- Put his "hat" on by pulling the floss through the slit at the top of his head
- next his "scarf" by winding the floss around the neck about 3 times
- then put his "mitten" on his 1st arm by winding the floss around the button once
- to put on his "boots" it is easiest to now turn him upside down, then pull the floss over his feet, threading the floss through the slots on the bottom of his feet
- it is best if you can hold his feet together a little while you bring the floss back up to his 2nd arm and put his other "mitten" on by winding the floss around the other button,
- then you go back up to the slit on the top of his head and finish by wrapping the extra floss back and forth from the slit to his neck.

Now Mr Flossman™ is ready to go "skiing", and he will need to "ski" around all of your teeth everyday!

It will be fun for both of you, just get in the habit!

